

Product Specifications



**AMERICAN  
BIOANALYTICAL**

Phone: (800)443-0600  
Fax: (508)655-2754  
Email: [info@americanbio.com](mailto:info@americanbio.com)

---

**Potassium Chloride**

**AB01652**

endotoxin tested

**Specifications**

CAS 7447-40-7  
grade Ultra Pure/USP  
formula weight 74.55  
formula KCl  
physical description white granular

	<b>Minimum</b>	<b>Maximum</b>
assay	99.0%	100.5%
identification	pass	
acidity or alkalinity	pass	
loss on drying		1.0%
iodide or bromide	pass	
aluminum		1µg/g
calcium and magnesium	pass	
sodium	pass	
heavy metals		0.001%
organic volatile impurities	pass	
endotoxin		2.5eu/g

**Notes and Applications**

DNase, RNase, Protease: none detected  
store at room temperature