



**AMERICAN  
BIOANALYTICAL**

Phone: **(800)443-0600**  
Fax: **(508)655-2754**  
Email: [info@americanbio.com](mailto:info@americanbio.com)

---

**Northern Hybridization Solution****AB10055**

DEPC treated

**Description**

a solution of 4X SSPE, 1% SDS, 5X Denhardt's Solution, 50% Formamide (deionized) and 5% Dextran Sulfate in DEPC treated water, 0.2µm filtered.

**Notes and Applications**

see AB13161, AB01922, AB03075, AB00600, AB000420, AB02128 for specifications  
store at -20°C

Reference:

Williams, S., Slatko, B., and McCarrey, J., *Laboratory Investigations in Molecular Biology*  
Jones and Bartlett, Sudbury, MA (2006)